



**Your trusted global partner in life-saving research**

For over 50 years, Allentown has provided critical solutions to the biomedical research community. Our high-quality equipment and advanced vivarium systems continue to define and exceed industry standards. Like everything we do, our unparalleled customer service is stamped with the Allentown hallmarks of integrity, dedication and care.


# Wi-Com Cloud™

**MONITORING OF CLOUD-ENABLED ECOFLO™ BLOWERS**



OVERVIEW

Providing critical solutions with integrity and care.



View Details	Blower Name	Blower Address	Firmware/App Version	Last Sync Date	ACH	Polarity
View Details	Blower A-IZ on desk	10.1.29.175	05.00 / V1.5	10/13/2022 09:21:52 AM	51	+
View Details		10.1.31.92	07.02 / V1.5	10/13/2022 09:08:02 AM	51	+
View Details		10.1.31.136	07.02 / V1.5	10/13/2022 09:07:32 AM	49	+
View Details	ADVProd1	192.168.14.107	05.01 / V1.5	10/13/2022 08:47:29 AM	51	+
View Details	4GB-OldTSLive	192.168.15.142	05.01 / V1.4	10/04/2022 12:24:55 PM	50	+
View Details	IzzyLab-EFS	10.1.29.35	05.00 / V1.4	10/03/2022 10:29:10 AM	61	+
View Details	Blower-B RP	10.1.2.207	05.12 / V1.4	08/18/2022 08:56:40 AM	50	+
View Details	BLOnRPDeskBlower	192.168.17.174	05.00 / V1.4	08/17/2022 02:01:22 PM	65535	+
View Details		10.1.28.215		08/11/2022 10:52:02 AM	65535	+
View Details	Winmate3 BL	192.168.17.186	05.12 / V1.4	07/14/2022 02:49:43 PM	50	+

## Digital Application for Monitoring Cloud-Enabled EcoFlo Blowers

The Wi-Com Cloud application allows you to monitor cloud-enabled EcoFlo blowers remotely or in person, in real time, by retrieving data via the app. With Wi-Com Cloud, data from multiple blowers across your facility will always be at your fingertips, ensuring that your facility is consistently operating smoothly.



**Continuous Monitoring** – Wi-Com Cloud allows you to monitor your blower's humidity level, airflow, rack exhaust pressure, filter status, and alarm status. Essentially, monitoring the full range of parameters and conditions that are necessary for peak IVC blower operation.

**User-Friendly Interface** – The Wi-Com Cloud user interface is intuitive and easy to use. No confusing or hard-to-understand processes, just simplicity.

**Alarm Notifications** – Wi-Com Cloud will help you stay in the know by sending you email or text notifications if there is a blower condition that requires attention. So, you'll always be prepared to quickly and effectively address issues before facility efficiency is impacted.

OVERVIEW

E...	Building	Room	Android Serial	IP Address	Blower Name	Customer
Edit	AlnBG-01	AFT Lab	200981D4D5AE5CF5	192.168.14.107	ADVProd1	Allentown
Edit	AlnBG-01	AFT Lab	170719D4F21FE5CF	192.168.17.174	BLOnRPDeskBlower	Allentown
Edit	AlnBG-01	AFT Lab	180621D4D5ACFFC8	10.1.29.51	Blower-D	Allentown
Edit	AlnBG-01	AFT Lab	182419D4F21FE5CF	192.168.1.231	IzzyHomeBL	Allentown
Edit	AlnBG-01	AFT Lab	081621D4D5ACFFC8	10.1.29.35	IzzyLab-EFS	Allentown
Edit	AlnBG-01	AFT Lab	1522W026603	192.168.17.186	Winmate3 BL	Allentown
Edit	AlnBG-01	AFT Lab	ABCDE12345	192.168.16.208	WinmateTestBL	Allentown
Edit	AlnBG-01	IT Lab	041779D4EE803CA3	192.168.15.142	4GB-OldTSLive	Allentown
Edit	AlnBG-01	IT Lab	0C2321D4D5ACFFC8	10.1.29.175	Blower-A-IZ on desk	Allentown
Edit	AlnBG-01	IT Lab	0B1E21D4D5ACFFC8	10.1.2.207	Blower-B RP	Allentown
Edit	AlnBG-01	IT Lab	090A21D4D5ACFFC8	10.1.28.252	Tower-S	Allentown



Features

- Android OS
- Wireless security WPA / WPA2 / WPA3-Enterprise
- LAN RJ45 connector
- Both static and DHCP IP addressing are available
- Enhanced system architecture
- Alarm notifications (email or email to text)
- Built-in antenna versus external
- Web versus client-based application
- Data collection and storage
- Reporting

Dashboard User Manage Facility Rack Blowers Reports

74 F 40% 51

07.02 / V1.5 10/13/2022 9:08:02 AM

Export to Excel

Blower Name	Blower Address	Sync Date	ACH	Polarity	Temp (S)	Humidity (S)	Temp (E)
No items to display							